

Zeiss AxioImager.Z2 Brightfield Protocol

1) System Startup


- ☞ Please note our sign-up policy. You must inform the facility at least 24 hours beforehand if you can't come; otherwise, you will receive a charge for unused time. The facility will allow for extenuating circumstances (cells dying, sick day, etc.) if you inform us in a timely fashion.
- ☞ Follow each step of the startup poster and wait for the microscope software to fully load.

2) Lens Cleaning

- ☞ Please **clean all of the lenses** (used and unused) **before and after** your session. Refer to the lens cleaning poster if you need any help recalling the rules and steps.


3) Microscope Control

- ☞ Take note of the microscope light path:
 - Light source
 - Direction through microscope
 - Condenser assembly
 - Slide
 - Lens
 - Camera/ Eyepiece
 - Note the locking positions of the sliders and remember to be delicate! The slider pulls glass.



- ☞ Touch pad control
 - Within the “Home” menu, select “Microscope” and then “Turret” (**Stay within this window for the entire session and ignore other tabs**)
 - Objective Menu: always begin focusing using the 10x lens.
 - (a) **Take care to properly set your focal position with the lower power lens before going to a higher power lens.** Proper focusing from low power to high power will help prevent the microscope from breaking your slide.

- Reflector Menu: select your first filter channel
- To open the shutter, turn the TL illumination on the far right of the screen “On”
 - (a) Remember to **close the shutter** as often as you open it to prevent rapid bleaching of your sample.
 - (b) Note: The active option is always white.



☞ Always begin imaging with the 10x objective and take care for the following issues:

- **Select and inspect each slide**
 - If it is dirty, gently clean with a Kim Wipes and/or cotton swab. You should do this with all your slides before you come.
- **Load the slide securely into the stage clip.**
 - Focusing will become difficult when the slide is uneven.

☞ Joystick Control

- Use the joystick to find your region of interest
- Use the F1 button in the upper right to toggle between course and fine x and y control



☞ Focus Control

- Course and fine focus control are located on both sides of the microscope.
- If you are focusing on the right side of the microscope, the fine focus is located **INSIDE** the coarse focus. You can adjust by using the finger positions within the coarse.
- Focus your sample and check all the channels before imaging.



4) Software Control


☞ Open the Zen Software

- Select the “Zen Pro” imaging option
- Click “Skip calibration”

☞ Upon Zen’s completed startup, please note the 3 divisions of the interface:

- Left: acquisition control
- Middle:
 - Top: image preview and acquired image


- Bottom: basic measurements and manipulations
- Right: image gallery



5) Your first image


- ☞ Select the H condenser filter and readjust light intensity using the LED dial
- ☞ Gently push both sliders all the way in to send the light to the eyepiece.
- ☞ Use the touchpad to:
 - **Select the 10x lens** in the objective menu
 - **Select the DIC/TL filter channel and open the TL illumination**
- ☞ Use the stage movement wand and focus knobs to **focus and center your specimen** on a region of interest.
- ☞ Setup Kohler illumination; you can refer to the poster for help. This step is critical and necessary for uniform field illumination in imaging and must be done every time you switch lenses.
- ☞ With one hand on focus and the other on the condenser iris for your desired level of brightness, contrast, and morphological detail.
 - Higher Lens NA → larger condenser aperture
 - Lower Lens NA → smaller condenser aperture

- ☞ Use the touchpad RL illumination option to **close the shutter**.
- ☞ **Pull the BOTTOM camera/eyepiece slider out** to send light to the camera.



1

2



3


4

Within the Zen software:


- ☞ Setup your experiment in the **Control Panel**:
- ☞ **AVOID THE LOCATE TAB!**
- ☞ Click the **“Acquisition”** tab.
- ☞ Select the **“Channels”** window.

- ☞ Enter the **“Live”** mode,
 - Click **“Set Exposure”** within the channels window.
- ☞ **Focus** on your sample using one of the following options:
 - Focus using the computer keyboard
 - Put your cursor over the image and hit **“Ctrl”** on the keyboard; use the mouse roller to fine tune the focus (to signify focus activation, a dialog box will appear).
 - Focus using the microscope fine focus
- ☞ Click the **“Acquisition Mode”** window to set the white balance.
 - Select either the auto option or picker option; find a white region of your sample and click using the picker option.

Before White Balance



After White Balance



- ☞ Click **“Stop”** (the same icon as **“Live”**)
- ☞ When all exposure times are set, Click **“Snap”** to acquire the image. Before saving your image **set your gamma = 1** using the display tab at the bottom of the display window.






Image Saving Options:

- ☞ Click the **“Save As”** button to save the image in your directory within the specified User Data folder.
 - Note: All data is stored in the D: drive under the “User Data” folder. Be sure to save to the “User Data” folder and NOT the “Users” folder.



- ☞ Optional: Auto-save function
 - Note: There is also an auto save function that must be enabled AND disabled before and after your session. You CANNOT export overlay images using this tool.




- Use the same parameters from the batch protocol to export your images correctly.

- Exporting options: be sure that you select the following image options to generate the same image you acquired.
 - (a) DO NOT select “Original Data” export
 - (b) Select “Apply Display Curve and Channel Color”



6) Higher magnification and switching slides

- ☞ If you switch to the higher magnification oil immersion lenses from a dry immersion lens, go back to the microscope and enter the objective menu on the touchpad.
 - Press the high power lens you want to use.
 - The lens will then move to the correct objective and the stage will drop.
 - The touchpad will prompt you to put immersion on the specimen.
 - At this point, apply **one drop of oil** and click “**Done**” to raise the stage back up.





- ☞ If you switch from an oil lenses to a dry lens, after the touchpad prompts you to remove the immersion, take the slide out, wipe off the oil with a cotton swab. Replace the slide into the stage clip and click **up arrow on the left** to raise the stage back up.

- ☞ When you need to switch slides,
 - Go back to the **10x** and refocus the new sample
 - If you are imaging a similar sample using the same objective lens, click the “**Load Position**” button located in the upper right hand corner of the touchpad. This will allow you to lower the stage, replace your slide with a new one, and restore the same focal plane.



7) Data Export

- Go to the **“Processing”** tab (located to the right of the **“Acquisition”** tab)
 - Click **“Batch”**
 - Under **“Batch Method”** select **“Image Export”**
 - In the middle of the screen under **“Batch Processing,”** click **“Add”** and then select all the images you want to export.
 - Highlight one image
 - Under parameters, make the following selections
 - File type: Tagged Image File.tiff
 - Compression: None
 - Unclick convert to 8 bit
 - Check **“Use Full Set of Dimensions”**
 - Unclick **“Create project folder”** unless you would like all images in separate folders
 - Exporting options: be sure that you select the following image options to generate the same image you acquired.
 - DO NOT select **“Original Data”** export
 - Select **“Apply Display Curve and Channel Color”**



- Underneath Batch Processing
 - Click **“Copy Parameters”**
 - Select all images. (Ctrl +Shift)
 - Click **“Paste Parameters”**
 - Click **“Run Selected”**

8) System Shutdown

- ☞ Back up all your data.
- ☞ Clean all the lenses.
- ☞ Check the microscope calendar to see when the next user has an appointment
- ☞ If the user comes within 2 hours, log off your account
 - Otherwise, follow shutdown poster steps.